



# MONOSCREEN<sup>®</sup> Ab ELISA

Instruction manual  
 BIOK192-NEO2\_NO\_(EN)\_V03  
 07/02/2025

## Monoscreen AbELISA *Neospora caninum* (SRS2)

Reference : BIO K 192

ELISA test for serodiagnosis of Bovine Neosporosis

Biwell, indirect test

For veterinary *in vitro* use only



Sample/Dilution	Bovine
Serum – Plasma* / 100X	✓
Milk / 4X	✓

\*Hereafter, we will refer to it as serum.

### Presentation

Product reference	BIO K 192/2	BIO K 192/5
Format	2 plates, strips of 16 wells	5 plates, strips of 16 wells
Reactions	96 tests	240 tests

### Kit composition

Provided material		BIO K 192/2	BIO K 192/5
Microplate	Microplates	2	5
Washing solution	Washing solution (20X)	1 x 100 mL	1 x 250 mL
Dilution solution	Colored dilution solution (5X)	1 x 50 mL	1 x 100 mL
TMB Solution	TMB Solution (1X)	1 x 25 mL	1 x 55 mL
Stop solution	Stopping solution (1X)	1 x 15 mL	1 x 30 mL
Conjugate	Conjugate (50X)	1 x 0,5 mL	1 x 1,4 mL
CTL POS	Positive control	1 x 0,5 mL	1 x 0,5 mL
CTL NEG	Negative control	1 x 0,5 mL	1 x 0,5 mL
Tracer	Tracer	1 x 0,5 mL	1 x 0,5 mL

### Revision history

Date	Version	Modifications
07/02/2025	V03	Layout and simplification of the entire manual

Note : minor typographical, grammar and formatting changes are not included in the revision history.


## A. Introduction

*Neospora caninum* is a protozoon that was originally described as a parasite in dogs, in which it causes myositis and encephalitis. Bovine neosporosis is now recognized as a major cause of spontaneous abortion in cattle. It is highly suspected on 20% of the farms with repeated abortions and a cow that is seropositive for *Neospora caninum* has a threefold greater risk of aborting than a cow that is *Neospora*-negative. *Neospora* is responsible for 21% of spontaneous abortions occurring in an individual animal. This percentage rises to 33% for the herd as a whole. Vertical transmission is the rule (at least 80% of the calves born to seropositive cows are infected). Serotesting before the calf's first colostrum intake will reveal prenatal infection.

## B. Test principle

The test uses 96-well microtitration plates sensitized by a purified *Neospora caninum* protein. The plate's odd columns (1, 3, 5, 7, 9, 11) contain the purified protein, whereas the even columns (2, 4, 6, 8, 10, 12) contain a control antigen. We thus have a genuine negative control. Using such a control reduces the number of false positives considerably.

The test blood sera, plasma or milks are diluted in the dilution solution. The plate is incubated and washed, then the conjugate, a peroxidase-labelled anti-bovine IgG1 monoclonal antibody, is added to the wells. The plate is then incubated a second time at  $21\pm 3^{\circ}\text{C}$  washed again and the chromogen tetramethylbenzidine (TMB) is added. This chromogen has the advantages of being more sensitive than the other peroxidase chromogens and not being carcinogenic. If specific anti-*Neospora caninum* immunoglobulins are present in the test sera or milks the conjugate remains bound to the microwell that contains the protozoon, and the enzyme catalyzes the transformation of the colorless chromogen into a pigmented compound. The intensity of the resulting blue colour is proportionate to the titre of specific antibody in the sample. The signal read off the negative control microwell is subtracted from that of the positive microwell sensitized by the protozoon protein.



## C. Material required but not provided

- Distilled/demineralized water.
- Dilution microplates (optional).
- Graduated mono or multichannel pipettes (2-20 $\mu\text{L}$ , 20-200 $\mu\text{L}$  and 10-1000 $\mu\text{L}$  range) and single-use tips.
- Microplate reader (450nm filter).
- Microplate washer (optional).
- Incubator at  $21\pm 3^{\circ}\text{C}$ .
- Standard laboratory equipment: graduated cylinder, tube rack, lid,...

## D. Warnings and precautions of use

- The reagents must be kept between  $+2$  and  $+8^{\circ}\text{C}$ .
- Unused strips must be stored with the desiccant in the hermetically sealed aluminum envelope.
- Do not use reagents beyond shelf-life date.
- Do not use reagents from other kits.
- Make sure to use distilled/demineralized water.

- The stopping solution contains 1M phosphoric acid. Handle it carefully.
- Used material must be disposed of in compliance with the legislation in force regarding environmental protection and biological waste management.
- Keep the TMB solution away from light.

## E. Preparation of the solutions

- The solutions are to be prepared extemporaneously.
- The washing solution must be diluted 20-fold in distilled/demineralized water. The cold solution crystallizes spontaneously. Bring the vial to  $21\pm 3^{\circ}\text{C}$  to make sure that all crystals have disappeared; mix the solution well and withdraw the necessary volume.
- The dilution solution is to be diluted 5 times in distilled/demineralized water. The dilution is colored in yellow. It is used for dilution of samples, kit controls (positive and negative controls, and tracer), and conjugate.
- The conjugate must be diluted 50-fold in the dilution solution.
- The stopping solution is ready to use.
- The TMB solution is ready to use. It must be perfectly colorless.
- The tracer is a reference sample that can be used to check the intra-laboratory reproducibility of the kit's batch. It must be diluted 100 times in the dilution solution.

## F. Preparation of the samples

- **Serum samples** and kit controls (positive and negative control, and tracer) are to be diluted **100 times** in the dilution solution and homogenized. Avoid using hemolyzed or coagulated samples.

*Recommended dilution:*

*10 $\mu\text{L}$  of sample + 990 $\mu\text{L}$  of dilution solution.*

- **Milk samples** are to be centrifuged 20 min at 4000g. Take up the middle layer of liquid taking care not to touch the underlying cell sediment. Samples are to be diluted **4 times** in the dilution solution and homogenized.

*Recommended dilution:*

*250  $\mu\text{L}$  of sample + 750 $\mu\text{L}$  of dilution solution.*

## G. Procedure

- Bring all the reagents to  $21\pm 3^{\circ}\text{C}$  before use.
  - Carefully read through the previous points.
1. Distribute the **diluted samples** and **diluted kit controls** at a rate of **100 $\mu\text{L}$  per well**. Cover and incubate the plate at  $21\pm 3^{\circ}\text{C}$  for **60 $\pm$ 5min**.
  2. Remove the content of the microplate. **Wash the microplate 3 times with 300 $\mu\text{L}$  of washing solution per well**. Avoid the formation of bubbles in the wells between each wash.
  3. Distribute the **diluted conjugate** at a rate of **100 $\mu\text{L}$  per well**. Cover with a lid and incubate at  $21\pm 3^{\circ}\text{C}$  for **60 $\pm$ 5min**.
  4. Remove the content of the microplate. **Wash the microplate 3 times with 300 $\mu\text{L}$  of washing solution per well**. Avoid the formation of bubbles in the wells between each wash.
  5. Distribute **100 $\mu\text{L}$  of TMB solution** per well. Incubate at  $21\pm 3^{\circ}\text{C}$  for **10 $\pm$ 1min** away from the light, without covering.
  6. Distribute the **stopping solution** at a rate of **50 $\mu\text{L}$  per well**. Color changes from blue to yellow.
  7. Record optical densities using a plate spectrophotometer with a **450nm filter** within **5 minutes** after adding the stopping solution.

## H. Validation of results

The test can only be **validated** if:

- The difference between the optical density (OD) readings of the odd and even wells of the positive control is greater than 0,800.

Positive control:  $OD_{\text{odd well}} - OD_{\text{even well}} > 0,800$

- The difference between the optical density (OD) readings of the odd and even wells of the negative control is less than 0,300.

Negative control:  $OD_{\text{odd well}} - OD_{\text{even well}} < 0,300$

## I. Interpretation of results

- Calculate for each sample its "delta OD" by subtracting the optical density from the even wells to the odd wells.
- Proceed the same way for positive and negative controls.

Divide each value obtained by the corresponding value obtained with the positive control and multiply this result by 100 to express it as a percentage.

$$\text{Val(ue)} = \frac{\text{Delta DO sample} * 100}{\text{Delta DO pos}}$$

By using the table below, determine the level of positive for the samples.

	Results	Status
Sample	$Val \leq 10\%$	Negative
	$10\% < Val \leq 15\%$	Doubtful
	$Val > 15\%$	Positive

These cut-off values apply to individual samples and diagnosis of abortive disease.

For optimal use, especially when testing mixtures of samples (bulk tank milk) or in epidemiological investigations, the cut-off of **6,5%** without a doubtful area is to be used.

Neg		Pos
$Val \leq$	6,5%	$< Val$

Get the interpretation of your results quickly and easily using **AnalysisScreen**, our free online platform, available on our website: <https://www.biox.com>.

# ANALYSISCREEN



AnalysisScreen™ is the new module for reading and interpreting all types of Monoscreen™ and Multiscreen™ ELISA plates. **AnalysisScreen™** is:


- Free
- Accessible online via our website: <https://www.biox.com>
- Updated in real time
- Compatible with all Bio-X Diagnostics plate designs
- Very easy to use



SCAN ME


Serum protocol

- 1 Dilution of samples 1/100  
Dilution of the kit controls (positive and negative control and tracer) 1/100



Milk protocol

- 1 Dilution of samples 1/4  
Dilution of the kit controls (positive and negative control and tracer) 1/100



Joint protocol

- 2 Add 100  $\mu$ L of diluted conjugate (1/50)



- 3 Add 100  $\mu$ L of TMB



- 4 Add 50  $\mu$ L of stop solution

- 5 Record the optical densities

450 nm



\* Notes do not replace the instructions for use of which they are a synthesis.