

## Monoscreen AbELISA Brucellosis

Reference : BIO K 140

ELISA Test for serodiagnostic of *Brucella abortus*

Monowell, indirect test

For veterinary *in vitro* only



Sample	Species	Individual analysis	Pool analysis *, possible up to
Serum/Plasma	Bovine	✓	10

\* This is done in accordance with the legislation in force in your country, the certifying body or the recommendations made by the NRL when they exist. Mixtures must be made volume to volume, i.e. by taking the same volume of each of the sera making up the mixture.

### Presentation

Product reference	BIO K 140/2	BIO K 140/5
Format	2 plates, strips of 8 wells	5 plates, strips of 8 wells
Reactions	192 tests	480 tests

### Composition of the kit

	Provided material	Type *	Code	BIO K 140/2	Code	BIO K 140/5
Microplates	Microplates	1	D01425	2	D01425	5
Washing solution (20X)	Washing solution (20X)	A	D00695	1 x 100 mL	D00696	1 x 250 mL
Dilution solution (1X)	Colored dilution solution (1X)	A	D01511	3 x 125 mL	D01555	3 x 250 mL
TMB solution (1X)	Single component TMB (1X)	A	D01585	1 x 30 mL	D01557	1 x 60 mL
Stop solution (1X)	Stop solution (1X)	A	D00680	1 x 30 mL	D01556	1 x 60 mL
Conjugate (50X)	Conjugate (50X)	1	D01596	1 x 0,6 mL	D01563	1 x 1,5 mL
CTL POS	Positive Control	a	D01396	1 x 0,5 mL	D01396	1 x 0,5 mL
CTL NEG	Negative Control	a	D01564	1 x 0,5 mL	D01564	1 x 0,5 mL

\*: (1): dependent on kit and batch : (a): dependent on kit / (A): substitutable with components A / (B): substitutable with components B.

### Revision history

Date	Version	Modifications
10/04/2025	V01	Creation
12/02/2026	V02	Addition of 2 plate packaging. Adjustment of component volume. Adjustment of validation results data.
23/04/2026	V03	Modification of kit name and introduction. Addition of table of symbols.

Note : minor typographical, grammar and formatting changes are not included in the revision history.

## A. Introduction

*Brucella* is the causative agent of brucellosis, an infectious and contagious disease in animals that can be transmitted to humans. The genus *Brucella* comprises ten species classified according to their pathogenicity and their hosts; seven species can be isolated from terrestrial mammals: *B. abortus*, *B. melitensis*, *B. suis*, *B. canis*, *B. ovnis*, *B. neotomae*, and *B. microti*.

The main animal reservoirs for *Brucella* are domestic cattle (*B. abortus*), sheep and goats (*B. melitensis*) and pigs (*B. suis*). In animals, brucellosis causes abortions, reduced fertility and a loss of milk production, which can result in significant economic losses. The disease is transmitted mainly through direct contact with infected animals or by consuming unpasteurized dairy products from infected animals. The bacteria can survive for several months outside the animal's body, in the external environment, particularly in cold and damp conditions. These bacteria in the environment remain a source of infection for other animals, which become infected through close contact (via the respiratory tract, conjunctiva or even by ingestion). Diagnosing brucellosis can be difficult as the symptoms are often similar to those of other diseases. Diagnostic tests are required to diagnose brucellosis.

## B. Test principle

This ELISA test has been developed and validated for the detection of *Brucella abortus* in cattle.

96-well microplates were sensitized with a synthetic *Brucella* antigen. Sera and controls were diluted in the dilution solution. After 60 minutes of incubation and a washing step, the operator adds the conjugate, peroxidase-coupled protein G. After a second 60-minute incubation and a second wash, the chromogen tetramethylbenzidine (TMB) is added. This chromogen has the dual advantage of being more sensitive than other peroxidase chromogens, and of not being carcinogenic.

If *Brucella*-specific immunoglobulins are present in the serum, the conjugate remains bound to the cup containing the *Brucella* antigen, and the enzyme catalyzes the transformation of the colorless chromogen into a blue product. The intensity of staining is proportional to the level of specific antibodies present in the sample.

## C. Material required but not provided

- Distilled/demineralized water.
- Dilution microplates.
- Graduated mono or multichannel pipettes (2-20  $\mu\text{L}$ , 20-200  $\mu\text{L}$  and 10-1000  $\mu\text{L}$  range) and single-use tips.
- Microplate washer (optional).
- Microplate reader (450 nm filter).
- Incubator at  $21\pm 3^\circ\text{C}$ .
- Standard laboratory equipment: graduated cylinder, tube rack, lid,...
- Dilution microplate.

### Additional kit

- **Tracer *Brucella* (Ref.: BDE K 140):** Internal reference material for Brucellosis serology by ELISA.

## D. Warnings and precautions of use

- The reagents must be kept between  $+2$  and  $+8^\circ\text{C}$ .
- Unused strips must be stored with the desiccant in the hermetically sealed aluminum envelope.
- Do not use reagents beyond shelf-life date.
- Make sure to use distilled/demineralized water.
- The stopping solution contains 1M phosphoric acid. Handle it carefully.
- Used material must be disposed of in compliance with the legislation in force regarding environmental protection and biological waste management.
- Keep the TMB solution away from light.

## E. Preparation of the solutions

- The solutions are to be prepared extemporaneously.
- The washing solution must be diluted 20-fold in distilled/demineralized water. The cold solution crystallizes spontaneously. Bring the vial to  $21\pm 3^\circ\text{C}$  to make sure that all crystals have disappeared; mix the solution well and withdraw the necessary volume.
- The dilution solution is ready to use. The dilution solution is colored yellow. It is used for the dilution of samples, kit controls, tracer and conjugate.
- The conjugate must be diluted 50-fold in the dilution solution.
- The stopping solution is ready to use.
- The TMB solution is ready to use. It must be perfectly colorless.

## F. Sample preparation

- Sera samples and kit controls (positive and negative) are to be diluted 1:100 in the dilution solution and homogenized. Avoid using hemolyzed or coagulated samples.

*Recommended two-step dilution:*

- 1) 10  $\mu\text{L}$  sample + 90  $\mu\text{L}$  dilution solution in dilution microplate.
- 2) 10  $\mu\text{L}$  of first step + 90  $\mu\text{L}$  of dilution solution in test plate.

## G. Procedure

- Bring all the reagents to  $21\pm 3^\circ\text{C}$  before use.
  - Carefully read through the previous points.
1. Distribute the **diluted serum samples** and the **diluted kit controls** at **100  $\mu\text{L}$  per well**. Cover and incubate plate at  **$21 \pm 3^\circ\text{C}$  for  $60 \pm 5$  min**.
  2. Remove the content of the microplate. **Wash the microplate 3 times** with **300  $\mu\text{L}$  of washing solution per well**. Avoid the formation of bubbles in the wells and the desiccation of the microplate between each wash.
  3. Add **100  $\mu\text{L}$  of diluted conjugate** per well. Cover and incubate plate at  **$21 \pm 3^\circ\text{C}$  for  $60 \pm 5$  min**.
  4. Remove the content of the microplate. **Wash the microplate 3 times** with **300  $\mu\text{L}$  of washing solution per well**. Avoid the formation of bubbles in the wells and the desiccation of the microplate between each wash.
  5. Distribute **100  $\mu\text{L}$  of TMB solution** per well. Incubate at  **$21 \pm 3^\circ\text{C}$  for  $10 \pm 1$  min** away from the light, without covering.
  6. Distribute the **stopping solution** at a rate of **100  $\mu\text{L}$  per well**. Color changes from blue to yellow.
  7. Record optical densities using a plate spectrophotometer with a **450 nm filter** within **5 minutes** of adding the stopping solution.

## H. Validation of results

The test can only be **validated** if:

- The difference between positive and negative control optical density readings is greater than 0,800.

$$\text{OD}_{\text{positive control}} - \text{OD}_{\text{negative control}} > 0,800$$

- The optical density of the negative control is less than 0,300.

## I. Interpretation of results

Calculate for each sample its coefficient (S/P %) using the following formula:

$$S/P (\%) = \frac{OD \text{ sample} - OD \text{ negative control}}{OD \text{ positive control} - OD \text{ negative control}} * 100$$

	Results	Status
Individual sample	%S/P < 40%	Negative
	%S/P ≥ 40%	Positive
Pool of 10	%S/P < 15%	Negative
	%S/P ≥ 15%	Positive

Get the interpretation of your results quickly and easily using AnalysiScreen, our free online platform, available on our website : <https://www.biox.com>



AnalysiScreen™ is the new module for reading and interpreting all types of Monoscreen™ and Multiscreen™ ELISA plates.

AnalysiScreen™ is :

- Free
- Accessible online via our website: <https://www.biox.com>
- Updated in real time
- Compatible with all Bio-X Diagnostics plate designs
- Very easy to use




SCAN ME

## Symbols

Symbol	Meaning
	Catalog number
	Manufacturer
	Temperature limitation
	Use by
	Batch code
	Consult Instructions for Use
	Contain sufficient for "n" tests
	Keep away from light
	Keep dry
	Corrosive substance
	Hazardous/irritating product

1 Distribute 100  $\mu$ L of diluted samples (1/100) and diluted kit controls (positive and negative control) (1/100)



2 Add 100  $\mu$ L of diluted conjugate (1/50)




3 Add 100  $\mu$ L of TMB



4 Add 100  $\mu$ L of stopping solution

5 Record optical densities



\* Notes do not replace the instructions of use of which they are a summary.